Pankreatische Amylase ELISA Kit

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Pancreatic amylase ELISA Kit

*For the in vitro determination of human pancreatic amylase in
stool*

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| | |
|---|------------|
| Table of content | Seite/Page |
| | 2 |
| 1. INTENDED USE | 13 |
| 2. SUMMARY AND EXPLANATION OF THE TEST | 13 |
| 3. PRINCIPLE OF THE TEST | 13 |
| 4. MATERIAL SUPPLIED | 14 |
| 5. MATERIAL REQUIRED BUT NOT SUPPLIED | 14 |
| 6. PREPARATION AND STORAGE OF REAGENTS | 15 |
| 7. PRECAUTIONS | 15 |
| 8. SPECIMEN COLLECTION AND PREPARATION | 16 |
| FAECES | 16 |
| 9. ASSAY PROCEDURE | 16 |
| PROCEDURAL NOTES | 16 |
| TEST PROCEDURE | 17 |
| 10. RESULTS | 18 |
| TYPICAL CALIBRATION CURVE | 18 |
| 11. LIMITATIONS | 19 |
| 12. QUALITY CONTROL | 19 |
| EXPECTED VALUES | 19 |
| 13. PERFORMANCE CHARACTERISTICS | 20 |
| SENSITIVITY | 20 |
| 14. REFERENCES | 20 |
| 15. GENERAL NOTES ON THE TEST AND TEST PROCEDURE | 21 |

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Lagerung/Storage: 2 - 8 °C

1. INTENDED USE

The *Immundiagnostik* Assay is intended for the quantitative determination of **Pancreatic Amylase** in stool. For in vitro diagnostic use only.

2. SUMMARY AND EXPLANATION OF THE TEST

Like lipase and elastase the **Pancreatic Amylase** is synthesized in the pancreas. In case of a chronic pancreatitis, **Pancreatic Amylase** is decreased. Some forms of disease (alcoholism, accidental trauma,) are associated with a pancreatic insufficiency. In the routine laboratory **Pancreatic Amylase** has proved itself to be a genuine alternative to Elastase-I in stool diagnosis.

Indication

- Chronic pancreatitis
- Exocrine pancreas insufficiency

3. PRINCIPLE OF THE TEST

In a first incubation step, the Pancreatic Amylase in the samples is bound to monoclonal mouse antibodies (in excess), which are immobilized to the surface of the microtitre wells. To remove all unbound substances, a washing step is carried out. In a second incubation step, an anti-Pancreatic Amylase (POD-monoclonal antibody) antibody is added. After another washing step, to remove all unbound substances, the substrate, tetramethylbenzidine (TMB) is added. An acidic stop solution is then added to stop the reaction. The colour converts from blue to yellow. The intensity of the yellow color is directly proportional to the concentration of Pancreatic Amylase in the sample. A dose response curve of the absorbance unit (optical density, OD) vs. concentration is generated, using results obtained from the calibrators. Pancreatic Amylase, present in the patient samples, is determined directly from this curve.

The combination of two specific antibodies in the Pancreatic Amylase ELISA drastically reduces the possibility of wrong-negatives results and offers a secure diagnostic system to the user.

4. MATERIAL SUPPLIED

| Catalogue No | Kit Components | Quantity |
|--------------|---|----------------|
| K 6410MTP | one holder with precoated strips | 12 x 8 |
| K 6410WB | ELISA wash buffer concentrate 10x | 2 x 100 ml |
| K 6410K | Konjugate (Maus anti-Pancreatic Amylase, peroxidase- labeled) | 1 x 30 μ l |
| K 6410ST | Calibrators, freeze dried (0; 440; 1750; 7000; 28000 mU/l) | 5 vials |
| K 6410 KO | Control, freeze dried | 1 vial |
| K 6410TMB | TMB substrate (Tetramethylbenzidine) | 1 x 15 ml |
| K 6410AC | ELISA stop solution, ready to use | 1 x 7 ml |

On demand we will send you 3 calibrator sets free of charge. Any further calibrator sets will be charged.

5. MATERIAL REQUIRED BUT NOT SUPPLIED

- Deionized water
- Precision pipettors calibrated to deliver 50-100 μ l
- A multi-channel dispenser or repeating dispenser
- Centrifuge capable of 3000 x g
- Vortex-Mixer
- Microplate reader 450 nm

6. PREPARATION AND STORAGE OF REAGENTS

- To run the assay more than one time, please make sure that the reagents are carefully stored as mentioned. Prepare just the appropriate amount necessary for the assay.
- The **ELISA wash buffer concentrate** should be diluted with Aqua dest. **1:10** before use (add 900 ml A. dest. to 100 ml concentrate). Crystals could occur due to high salt concentration. The crystals have to be resuspended **before dilution of the buffer solutions** using a water bath (37°C). The buffer concentrates are stable at 2-8°C until the expiry date stated on the label. Diluted solutions could be stored at 2-8°C for 1 month.
- The **standards** and the **control** have to be reconstituted with **250 µl aqua dest.** Allow the vials to stand for **10 minutes** and then mix thoroughly by gentle inversion to insure complete reconstitution.
- The **POD antibody** has to be diluted **1:1000** in wash buffer (10 µl POD antibody and 10 ml wash buffer). The antibody is stable at 2 -4 °C until expiry date given on the label. **Diluted antibody solution is not stable and could not be stored.**
- All other test reagents are ready for use. The test reagents are stable up to the date of expiry (see label of test package) when stored at 4°C.

7. PRECAUTIONS

- For in vitro diagnostic use only.
- The calibrators and controls contain human source material which was tested and found to be non-reactive to HBsAg, anti-HIV-1/2, and anti-HCV. Since no method can offer complete assurance that hepatitis B virus, HIV-1/2, HVC or other infectious agents are absent, these reagents should be handled as if potentially infectious.
- Stop Solution consists of Sulfuric Acid. This is a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped out immediately with copious quantities of water. Do not breath vapor and avoid inhalation.
- Reagents should not be used beyond the expiration date shown on kit label.

8. SPECIMEN COLLECTION AND PREPARATION

Faeces

Give about **100 mg** of the sample (please note the real weight for the calculation) to **5 ml** of the wash buffer and mix. Centrifuge the diluted sample for 10 min at 3000 rpm. 1 ml supernatant is given into an eppendorf tube and centrifuged once more at 13.000 rpm for 5 min.

Afterwards the supernatant is diluted **1:40** in wash buffer (25 µl supernatant is added to 975 µl wash buffer). **100 µl** of this solution is used for the assay.

We recommend to weight out the stool samples for each run new. Supernatant is not stable and can't be stored.

Stool samples can be stored at -20°C for 4 weeks. Avoid thawing and freezing cycles.

Immundiagnostik recommends the use of the sample tubes from Roche Diagnostics for sample preparation.

9. ASSAY PROCEDURE

Procedural notes

- Do not interchange different lot numbers of any kit component within the same assay.
- The quality control guidelines should be observed.
- Incubation time, incubation temperature and pipetting volumes of the different components have been defined by the producer. Any variations of the test procedure, without consulting the producer with the producer, may influence the test results. Immundiagnostik AG can therefore not be held responsible for any damage.
- Conduct the assay according to the manual delivered with the kit.

Test procedure

Wash the precoated microtiter plate 5 x with 250 µl ELISA wash buffer. Carry out the tests in duplicate.


1. Add **100 µl** standard solutions and prediluted patient samples
2. Incubate for **1 hour**, shaking on a horizontal mixer, at room temperature.
3. Decant the content of the plate and wash the wells **5 x with 250µl** ELISA wash buffer.
4. Add **100 µl** prediluted anti-Pancreatic Amylase antibody.
5. Incubate for **1 hour**, shaking on a horizontal mixer, at room temperature.
6. Decant the content of the plate and wash the wells **5 x with 250µl** ELISA wash buffer.
7. Add **100 µl** TMB substrate solution.
8. Incubate for **5-15 minutes** at room temperature.
9. Add **50 µl** stop solution and mix shortly.
10. Determine absorption with an ELISA reader at **450 nm** against 620 nm as reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the measurement range of the photometer, absorption must be measured immediately at 405 nm against 620 nm as reference.

10. RESULTS

A calibration curve is constructed from the standards. Commercially available software can be used as well as graph paper. Results of the samples are read from this calibration curve.

THE CALIBRATION CURVE IS NOT LINEAR, therefore a spline- or 4PL algorithm is recommended.

Typical calibration curve

| | | | | | |
|----------------------|------|-------|-------|-------|-------|
| Konzentration [mU/l] | 0 | 440 | 1750 | 7000 | 28000 |
| OD Mittelwert | 0.05 | 0.195 | 0.651 | 1.812 | 2.544 |

These data are for demonstration only and cannot be used instead of data obtained from the actual assay

Faeces

For the Pancreatic Amylase concentration of faeces samples, calculate as described in the following example:

weight: 80 mg (1ml Stool = 1g) = 0,08 ml

dilution step 1: 5ml / 0,08ml = 62,5

dilution step 2: 40

dilution factor: 62,5 x 40 = 2500

Multiplicate the result with **2500** to get the real concentration. The dilution factor depends on the weight of the faeces.

11. LIMITATIONS

Samples with Pancreatic Amylase levels greater than the highest calibrator, should be diluted and re-assayed.

12. QUALITY CONTROL

Immundiagnostik recommends commercial control samples for internal quality control.

Control samples or serum/stool pools should be analyzed with each run of calibrators and patient samples. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. In assays in which one or more of the quality control sample values lie outside the acceptable limits, the results for the patient sample may not be valid.

Expected values

| | |
|-----------------------------|-----------------|
| Pancreatic Amylase (stool): | > 1000 U/l |
| Grey area | 1000 – 1500 U/l |

13. PERFORMANCE CHARACTERISTICS

Sensitivity

The detection limit was defined as $B_0 + 2SD$.

n=20

| Sample | Pancreatic Amylase Mean value [OD] | Standard variation | Detection limit [mU/l] |
|--------|------------------------------------|--------------------|------------------------|
| 1 | 0.065 | 0.030 | 237 |

14. REFERENCES

1. Bishop et al.: 1996; Pancreas, 13 (3), 226
2. Katschinsky et al.: 1997; Pancreas, 15 (2), 1991

15. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- The test components which are made of human serum are tested for Australia antigen and HIV and found to be negative. However, since no test method can offer complete assurance that infectious agents are absent, these reagents should be handled as recommended for any potentially infectious human serum or blood specimen. The normal precautions for laboratory working should be observed.
- Reagents of the test package contain sodium azide as a bactericide. Contact with skin or mucous membranes has to be avoided.
- All reagents in the test package are to be used for in-vitro diagnostics only.
- The reagents should not be used after the date of expiry (see label on the test package).
- Single components with different lot numbers should not be mixed or exchanged.
- The guidelines for medical laboratories should be observed.
- Incubation time, incubation temperature and pipetting volumes of the different components are defined by the producer. Any variations of the test procedure, that are not coordinated with the producer, may influence the results of the test. Immundiagnostik can therefore not be held reliable for any damage resulting from this.