Osteoprotegerin (human) ELISA Kit

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- Jung K et al. (2001) Osteoprotegerin in serum as a novel marker of bone metastatic spread in prostate cancer. *Clin Chem* 47(11):2061-63
- Szulc P et al. (2001) Osteoprotegerin serum levels in men: correlation with age, estrogen, and testosterone status. *J Clin Endocrinol Metab* 86(7):3162-65
- Kudlacek S et al. (2003), Serum levels of osteoprotegerin increase with age in a healthy adult population. *Bone* 32: 681-686

Osteoprotegerin allgemein

- Hofbauer LC (1999), Osteoprotegerin ligand and osteoprotegerin: novel implications for osteoclast biology and bone metabolism. *European Journal of Endocrinology* 141: 195-210
- Aubin JE et al. (2000) Osteoprotegerin and its Ligand: A New Paradigm for Regulation of Osteoclastogenesis and Bone Resorption. *Medscape Women Health* 5(2)
- Simonet WS et al. (1997) Osteoprotegerin, a novel secreted protein involved in the regulation of bone density. *Cell* 89: 309-319
- Bucay N et al. (1998) Osteoprotegerin-deficient mice develop early onset osteoporosis and arterial calcification. *Genes and Development* 12: 1260- 1268
- Bekker PJ et al. (1999) Osteoprotegerin (OPG) has Potent and Sustained Anti-Resorptive Activity in Postmenopausal Women. 21st Annual Meeting of the ASBMR Abstract No. 1190
- Makhluף HA et al. (2000) Age-Related Decline in Osteoprotegerin Expression by Human Bone Marrow Cells Cultured in Three-Dimensional Collagen Sponges. *Biochemical and Biophysical Research Communications* 268: 669-672
- Honore P et al. (2000), Osteoprotegerin blocks bone cancer-induced skeletal destruction, skeletal pain and pain-related neurochemical reorganization of the spinal cord. *Nature Medicine* 6(5): 521-528

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Hersteller

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
Osteoprotegerin (human) ELISA Kit

For the in vitro determination of OPG (human) in serum,
plasma, urine and cell culture supernatant

Gültig ab / Valid from 29.02.2008

REF KB 1011

IVD

Immundiagnostik AG

Content

1. INTENDED USE	15
2. SUMMARY AND EXPLANATION OF THE TEST	15
3. PRINCIPLE OF THE TEST	15
4. MATERIAL SUPPLIED	16
5. MATERIAL REQUIRED BUT NOT SUPPLIED	17
6. PREPARATION AND STORAGE OF REAGENTS	17
7. PRECAUTIONS	18
8. SPECIMEN COLLECTION AND PREPARATION	18
9. ASSAY PROCEDURE	18
Procedural notes	18
Test procedure	19
10. RESULTS	20
11. LIMITATIONS	20
12. QUALITY CONTROL	20
Expected values	20
13. PERFORMANCE CHARACTERISTICS	22
Precision and reproducibility	22
Sensitivity	22
14. REFERENCES	23
15. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	25

1. INTENDED USE

The Immundiagnostik assay is a sandwich ELISA intended for the quantitative determination of Osteoprotegerin in serum, plasma, urine and cell culture supernatants. It is for in vitro diagnostic use only.

2. SUMMARY AND EXPLANATION OF THE TEST

Osteoprotegerin (OPG) or **Osteoclastogenesis inhibitory factor (OCIF)** is a dimeric glycoprotein of the TNF receptor family with a molecular weight of 60 kD resp. 120 kD which shows an inhibitory effect on osteoclasts and osteoclast precursor cells.

Osteoprotegerin is a soluble "decoy"-receptor and is produced in different tissues, e.g. bone, skin, liver, stomach, intestine and lung. As a so-called "decoy receptor" OPG inhibits the binding of RANK to RANKL (OPG-L, osteoclast differentiation factor, ODF) and thus inhibits the recruitment, proliferation and activation of osteoclasts.

OPG shows an inhibitory effect on osteoclasts. Osteoclast formation activity may be determined principally by the relative concentration of OPG-L/osteoclast differentiation factor (ODF) to OPG/OCIF in the bone marrow microenvironment. Alterations of this ratio may be the major cause of bone loss in many imbalances in bone metabolism such as osteoporosis, osteopetrosis, metastatic osteolytic lesions and rheumatic bone degradation.

Indication

- Postmenopausal and senile osteoporosis
- Glucocorticoid-induced osteoporosis
- Diseases with locally increased resorption activity
- Therapy monitoring after treatment with OPG
- Arthritis
- Oncology

3. PRINCIPLE OF THE TEST

This sandwich-type ELISA is an assay for the direct determination of OPG in serum, plasma and urine. In this assay two highly specific antibodies against OPG are used. The binding antibody is attached to the wells of the microtiterplate, the detection antibody is labeled with biotin.

In a first incubation step the samples and the biotinylated antibody against OPG react simultaneously with the pre-coated antibody on the microtiterplate. A sandwich-type complex is found consisting of the binding antibody on the plate, OPG and the biotinylated detec-

tion antibody. To remove all unspecific bound substances a washing step is carried out. In a second step streptavidin – peroxidase is added which reacts with the detection antibody. After another washing step the solid phase is incubated with the substrate, TMB. An acidic stopping solution is subsequently added. The blue colour changes to yellow. The intensity of the yellow colour is directly proportional to the concentration of OPG in the sample.

A dose - response curve of the absorbance units (at 450 nm) versus concentration is generated. OPG, present in the samples, is determined directly from this calibration curve.

4. MATERIAL SUPPLIED

Catalogue no.	Content	Kit Components	Quantity
KB1011MTP	PLATE	Pre-coated microtiterplate, 12 x 8 strips	96
KB1011WP	WASHBUF	ELISA washing buffer concentrate (20x)	50 ml
KB1011AP	ASYBUF	Assay buffer, ready-to-use	25 ml
KB1011ST	STD	Calibrators, ready-to-use (0; 0.37; 1.1; 3.3; 10; 30 pmol/l)	6 x 500 µl
KB1011ST	STDKONZ	Stock solution (500 pmol/l)	500 µl
KB1011KO	CTRL	Control, ready-to-use	500 µl
KB1011A2	AB	2 nd antibody (anti-OPG, biotinylated), ready-to-use	7 ml
KB1011K	CONJ	Conjugate, (streptavidin-HRP-labeled), ready-to-use	22 ml
KB1011TMB	SUB	TMB substrate (Tetramethylbenzidine), ready-to-use	22 ml
KB1011AC	STOP	ELISA stop solution, ready-to-use	7 ml

5. MATERIAL REQUIRED BUT NOT SUPPLIED

- 1,5 ml reaction vials (Eppendorf)
- Precision pipettes calibrated to deliver 10 -1000 µl and disposable tips.
- Centrifuge capable of 3000 x g
- ELISA reader
- Vortex-mixer
- Bidistilled or deionized water

6. PREPARATION AND STORAGE OF REAGENTS

- To run assay more than once, ensure that reagents are stored at conditions stated on the label. **Prepare only the appropriate amount necessary for each assay.** The kit can be used up to 4 times within the expiry date stated on the label.
- Reagents with a volume less than **100 µl** should be centrifuged before use to avoid loss of volume.
- The **ELISA WASHBUF** (wash buffer concentrate) should be diluted with aqua bidist. **1:20** before use (50 ml concentrate + 950 ml aqua bidist.), mix well. Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved at room temperature or at 37°C before dilution. The **WASHBUF** is stable at **2-8°C** until the expiry date stated on the label. Diluted **buffer solution** can be stored in a closed flask at **2-8°C for one month.**

▪ Cell culture

The **STDKONZ** (OPG calibrator concentrate, 500 pmol/l = 10 000 pg/ml) must be diluted **1:20** with medium (50 µl Calibrator + 950 µl Medium) to obtain **S1**. Allow the vial content to dissolve for 10 minutes and mix thoroughly by gentle inversion to insure complete reconstitution. Standard curve solutions must be prepared from S1 in **1:2** dilution steps by adding dilution medium as follows:

			S1	(25.0 pmol/l = 500 pg/ml)
200 µl S1	+	200 µl medium	= S2	(12.5 pmol/l = 250 pg/ml)
200 µl S2	+	200 µl medium	= S3	(6.25 pmol/l = 125 pg/ml)
200 µl S3	+	200 µl medium	= S4	(3.13 pmol/l = 62.6 pg/ml)

Cell culture medium is used as Standard 0.

7. PRECAUTIONS

- For in vitro diagnostic use only.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Stop solution is composed of sulfuric acid, which is a strong acid. Even diluted, it still must be handled with care. It can cause acid burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spills should be wiped out immediately with copious quantities of water.
- Reagents should not be used beyond the expiration date stated on kit label.

8. SPECIMEN COLLECTION AND PREPARATION

Serum, plasma and urine samples:

Serum, plasma and urine samples can be used **without any dilution**.

Serum must be centrifuged and aliquoted within 90 min after collection and stored at $-20\text{ }^{\circ}\text{C}$ until use.

Lipton et al. (2002) Serum Osteoprotegerin Levels in Healthy Controls and Cancer Patients. *Cancer Res.* 8:2306-2310

9. ASSAY PROCEDURE

Procedural notes

- Do not mix different lot numbers of any kit component.
- Guidelines for medical laboratories should be followed.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- The assay should always be performed according the enclosed manual.

Test procedure

<p>It is recommended to carry out the assay in duplicate and to cover the microtiter plate during the incubation steps.</p>
1. Add 100 µl ASYBUF (assay buffer) into each well..
2. Add 50 µl STD (standard), CTRL (control) or samples into each microtiter well in duplicate.
3. Add 50 µl AB (2 nd antibody) into each well. Mix gently.
4. Incubate over night (18 – 24 h) at 2-8 °C.
5. Aspirate and wash the wells 5 x with 300 µl ELISA wash buffer.
6. Add 200 µl CONJ (conjugate) into each well.
7. Incubate for 1 hour at room temperature (18 – 26 °C).
8. Aspirate and wash the wells 5 x with 300 µl ELISA wash buffer.
9. Add 200 µl SUB (TMB substrate solution) into each well.
10. Incubate for 15 – 20 min at room temperature in the dark.
11. Add 50 µl STOP (stop solution) into each well and mix shortly.
12. Determine absorption immediately with an ELISA reader at 450 nm against 620 nm (or 690 nm) as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at 405 nm against 620 nm as a reference.

10. RESULTS

A calibration curve is constructed from the standards. Commercially available software can be used as well as graph paper. Results of the samples are calculated from this calibration curve.

THE CALIBRATION CURVE IS NOT LINEAR, therefore a 4PL- algorithm or spline is recommended.

11. LIMITATIONS

Samples with OPG levels greater than the highest standard value should be further diluted with wash buffer and re-assayed.

12. QUALITY CONTROL

Immundiagnostik recommends to use control samples for internal quality control.

Control samples or serum pools should be analyzed with each run of calibrators and patient samples. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. In assays in which one or more of the quality control sample values are located outside the acceptable limits, the results for the patient sample may not be valid.

Expected values

Reference values:

Normal sera measured with our OPG-ELISA gave the values listed in the table (see next page):

Table (Kudlacek S et al., 2003):

OPG Reference data (n = 1134, 50.83 ± 51.47 pg/ml, median 36, 2-584;
n = 1134, 2.54 ± 2.57 pmol/L, median 1.8, 0.1-29.2)

Age	OPG (mean ± SD)		OPG (median, range)		n
	pg/ml	pmol/L	pg/ml	pmol/L	
Female population (n = 687, 51.9 ± 51.3 pg/ml, median 39; 2.60 ± 2.57 pmol/L, median 1.95)					
< 30	44.5 ± 21.2	2.23 ± 1.06	41.7, 87.7	2.09, 4.39	48
31-40	44.1 ± 1.25	2.21 ± 1.69	36, 260	1.8, 13	198
41-50	41.2 ± 25.0	2.06 ± 33.8	36, 204	1.8, 10.2	217
51-60	39.5 ± 22.3	1.98 ± 1.12	35.5, 160	1.78, 8	150
61-70	69.6 ± 62.9	3.48 ± 3.15	60.5, 342	3.03, 17.1	34
71-80	134.0 ± 70.0	6.70 ± 3.50	131, 244	6.55, 12.2	16
> 81	227.0 ± 100.0	11.35 ± 5.00	206, 479	10.3, 23.95	24
Male population (n = 447, 50 ± 51.7 pg/ml, median 34.8; 2.5 ± 2.59 pmol/L, median 1.74)					
< 30	40.7 ± 26.3	2.04 ± 1.32	32, 96	1.6, 4.8	19
31-40	41.4 ± 30.5	2.07 ± 1.53	33.5, 170	1.68, 8.5	72
41-50	36.3 ± 31.8	1.82 ± 1.59	36.2, 316	1.81, 15.8	116
51-60	36.9 ± 20.2	1.85 ± 1.01	34, 134	1.7, 6.7	176
61-70	41.0 ± 22.7	2.05 ± 1.14	29.8, 160	1.49, 8	32
71-80	119.0 ± 71.0	5.95 ± 3.55	38, 103	1.9, 5.15	12
> 81	226.0 ± 78.0	11.30 ± 3.90	130, 219	6.5, 10.95	20

Conversion factor: 1 pg/ml = 0.05 pmol/l

We recommend each laboratory to establish its own norm concentration range

13. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

The precision (intra-assay variation) of the Immundiagnostik OPG ELISA test was calculated from 16 replicate determinations on each of one samples.

Table1: **Intra-Assay** CV (n= 16)

Sample	OPG [pmol/l]	Intra-Assay CV [%]
1	4.5	10
2	19.5	8

The total precision (inter-assay variation) was calculated from data on 2 samples obtained in 16 different assays by three technicians on two different lots of reagents over a period of three months.

Table 2: **Inter-Assay** CV (n= 16)

Sample	OPG [pmol/l]	Inter-Assay CV [%]
1	4.5	< 10
2	19.5	< 10

Sensitivity

The detection limit (0 pg/ml + 3 SD): 0.14 pmol/l

14. REFERENCES

Publications based on Immundiagnostik's Osteoprotegerin ELISA:

- Simonini G et al. (2005) Osteoprotegerin serum levels in Kawasaki disease: an additional potential marker in predicting children with coronary artery involvement. *J Rheumatol* 32(11):2233-38
- Franck H et al. (2004) Evaluation of bone mineral density, hormones, biochemical markers of bone metabolism, and osteoprotegerin serum levels in patients with ankylosing spondylitis. *J Rheumatol* 31(11):2236-41
- Hofbauer LC et al. (2004) Effects of oral contraceptives on circulating osteoprotegerin and soluble RANK ligand serum levels in healthy young women. *Clin Endocrinol (Oxf)* 60(2):214-19
- Schoppet M et al. (2003) Increased osteoprotegerin serum levels in men with coronary artery disease. *J Clin Endocrinol Metab.* 88(3):1024-28
- Avbersek-Luznik I et al. (2002) Increased levels of osteoprotegerin in hemodialysis patients. *Clin Chem Lab Med* 40(10):1019-23
- Coen G et al. (2002) Serum osteoprotegerin and renal osteodystrophy. *Nephrol Dial Transplant* 17(2):233-38
- Jung K et al. (2002) Osteoprotegerin and receptor activator of nuclear factor-kappaB ligand (RANKL) in the serum of healthy adults. *Int J Biol Markers* 17(3):177-81
- Liegibel UM et al. (2002) Concerted action of androgens and mechanical strain shifts bone metabolism from high turnover into an osteoanabolic mode. *J Exp Med* 196(10):1387-92
- Lipton et al. (2002) Serum Osteoprotegerin Levels in Healthy Controls and Cancer Patients. *Cancer Res.* 8:2306-2310
- Sattler AM et al. (2002) Osteoprotegerin in patients with coronary artery disease. Abstract (Poster 002) of XIII Lipid Meeting, Leipzig
- Viereck V et al. (2002) Bisphosphonates pamidronate and zoledronic acid stimulate osteoprotegerin production by primary human osteoblasts. *Biochem Biophys Res Commun* 291(3):680-86
- Viereck V et al. (2002) Phytoestrogen genistein stimulates the production of osteoprotegerin by human trabecular osteoblasts. *J Cell Biochem* 84(4):725-35
- Burguera B et al. (2001) Leptin reduces ovariectomy-induced bone loss in rats. *Endocrinology* 142(8):3546-53
- Jung K et al. (2001) Osteoprotegerin in serum as a novel marker of bone metastatic spread in prostate cancer. *Clin Chem* 47(11):2061-63

- Szulc P et al. (2001) Osteoprotegerin serum levels in men: correlation with age, estrogen, and testosterone status. *J Clin Endocrinol Metab* 86(7):3162-65
- Kudlacek S et al. (2003), Serum levels of osteoprotegerin increase with age in a healthy adult population. *Bone* 32: 681-686

Osteoprotegerin general:

- Hofbauer LC (1999), Osteoprotegerin ligand and osteoprotegerin: novel implications for osteoclast biology and bone metabolism. *European Journal of Endocrinology* 141: 195-210
- Aubin JE et al. (2000) Osteoprotegerin and its Ligand: A New Paradigm for Regulation of Osteoclastogenesis and Bone Resorption. *Medscape Women Health* 5(2)
- Simonet WS et al. (1997) Osteoprotegerin, a novel secreted protein involved in the regulation of bone density. *Cell* 89: 309-319
- Bucay N et al. (1998) Osteoprotegerin-deficient mice develop early onset osteoporosis and arterial calcification. *Genes and Development* 12: 1260- 1268
- Bekker PJ et al. (1999) Osteoprotegerin (OPG) has Potent and Sustained Anti-Resorptive Activity in Postmenopausal Women. 21st Annual Meeting of the ASBMR Abstract No. 1190
- Makhluף HA et al. (2000) Age-Related Decline in Osteoprotegerin Expression by Human Bone Marrow Cells Cultured in Three-Dimensional Collagen Sponges. *Biochemical and Biophysical Research Communications* 268: 669-672
- Honore P et al. (2000), Osteoprotegerin blocks bone cancer-induced skeletal destruction, skeletal pain and pain-related neurochemical reorganization of the spinal cord. *Nature Medicine* 6(5): 521-528

15. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and distributed according to the IVD guidelines of 98/79/EC.
- All reagents in the kit package are for in vitro diagnostic use only.
- Guidelines for medical laboratories should be followed.
- Human materials used in kit components were tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Kit reagents contain sodium azide or thimerosal as bactericides. Sodium azide and thimerosal are toxic. Substrates for the enzymatic color reactions are toxic and carcinogenic. Avoid contact with skin or mucous membranes.
- Reagents should not be used beyond the expiration date shown on the kit label.
- Do not mix different lot numbers of any kit component.
- Incubation time, incubation temperature and pipetting volumes of the components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held responsible for any damage resulting from wrong use.
- Warranty claims and complaints in respect of deficiencies must be logged within 14 days after receipt of the product. The product should be send to Immundiagnostik AG along with a written complaint.

This Assay was developed by
Immundiagnostik, Bensheim and Biomedica, Wien.

Used Symbols:

Store at

Catalog Number

In Vitro Diagnostic Device

No. of tests

Manufacturer

Use by

Lot number

Immundiagnostik AG

Stubenwald-Allee 8a
D-64625 Bensheim

Tel.: +49(0) 62 51/70 19 00

Fax: +49(0) 62 51/84 94 30

info@immundiagnostik.com

www.immundiagnostik.com