Coenzym Q₁₀

HPLC

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*For the determination of Coenzyme Q₁₀ (Ubiquinone) in
EDTA-blood, serum and EDTA-plasma*

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1. INTENDED USE	15
2. SUMMARY AND EXPLANATION OF THE TEST	15
3. PRINCIPLE OF THE TEST	16
4. MATERIAL SUPPLIED	17
5. MATERIAL REQUIRED BUT NOT SUPPLIED	17
6. PREPARATION AND STORAGE OF REAGENTS	18
7. PRECAUTIONS	18
8. SPECIMEN COLLECTION AND PREPARATION	18
9. ASSAY PROCEDURE	19
PROCEDURAL NOTES	19
SAMPLE AND STANDARD PREPARATION	19
CHROMATOGRAPHIC CONDITIONS	20
10. TREATMENT OF THE COLUMN	20
11. RESULTS	20
CALCULATION	20
TYPICAL CHROMATOGRAM	21
12. LIMITATIONS	21
13. QUALITY CONTROL	21
EXPECTED VALUES	21
CONTROLS	22
14. PERFORMANCE CHARACTERISTICS	22
PRECISION AND REPRODUCIBILITY	22
LINEARITY	22
DETECTION LIMIT	22
15. DISPOSAL	22
16. TROUBLESHOOTING	23
17. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	24

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1. INTENDED USE

The *Immundiagnostik* Assay is intended for the quantitative determination of Coenzyme Q10 (Ubiquinone) in EDTA-blood, serum and EDTA-plasma. This Assay is designed for in vitro diagnostic use only.

2. SUMMARY AND EXPLANATION OF THE TEST

Ubiquinone was first isolated in the 50th from the group of Prof. Green (Wisconsin). The function was investigated by Prof. Mitchel, who received the Nobel price for his research on the oxydative phosphorylation pathway.

Ubiquinone is a coenzyme which is represented in every cell and the whole metabolism. It is built up by a chinonic ring and an isoprenic sidechain. In humans ubiquinone can be synthesized and taken up by nutrition.

Ubiquinone has two different pysiological functions.

- 1.) Component of the energy metabolism
- 2.) Radical scavanger

During the reduction of oxygen in the oxydative phosphorylation 3 Mol ATP are generated. In this reduction electrons are transferred from NADPH to oxygen via 6 different redox systems. Ubiquinone is the less abundant redox system in the membrane of the mitochondria. Because of the low amount it is the speed controlling redox-system in the energy metabolism. Normally the amount of ubiquinone is sufficient, but with growing age and exposure to sunlight it is reduced to 50 %.

Ubiquinone has a high amount of carbon doublebonds and therefore a higher potential of reduction than vitamin C or vitamin E. Thereby it is the first line of defense against free radicals. Therefore ubiquinone is an optimal stabilizer of the ion channels of the membranes.

Indications

- Determination of Ubiquinone status
- Cardiovascular disease
- Carcinogenesis
- Aging
- Burnout syndrome

3. PRINCIPLE OF THE TEST

The first step in determining Coenzyme Q10 includes sample preparation. Therefore the internal standard is added to the calibrator, controls and samples. Higher molecular substances are removed by precipitation and centrifugation. The Coenzyme Q10 in the supernatant is then extracted by an organic solvent, which is evaporated afterwards. The residue is resuspended in ethanol p.A.. 100 µl of the resuspension are injected into the HPLC system.

The separation via HPLC follows an isocratic method at 30 °C using a reversed phase column. One run lasts 15 minutes. The chromatograms are recorded by a UV detector. The quantification is performed with the delivered blood calibrator; the concentration is calculated via integration of the peak areas by the internal standard method.

The kits includes all reagents for preparation and separation of the samples, except the column.

Summary:

Besides many other parameters the advantage of HPLC method lies in the simultaneous handling of many analytes in a single test. The HPLC system enables even laboratories without experience in high performance liquid chromatography to use this technique for clinical routine determination in a quick and precise manner. Unlike immuno assays with up to six calibrators per test, a one-point calibration is sufficient to calibrate the test system. It is possible to automate the sample application and calculation of the results so that even higher numbers of samples can be handled nearly without control.

4. MATERIAL SUPPLIED

Catalogue No	Kit Components	Quantity
KC 1700lm	Mobile phase	1 x 1000 ml
KC 1700ka	Calibrator (lyophilized 0.5 ml)	5 vial
KC 1700is	Internal standard	1 x 110 ml
KC 1700rk	Reconstitution solution	1 x 5 ml
KC 1700vl	Dilution solution	1 x 85 ml
KC 1700ex	Extraction solution	1 x 220 ml
KC 1700et	Ethanol p.A.	1 x 20 ml

5. MATERIAL REQUIRED BUT NOT SUPPLIED

- 10 ml screw capped glas vials (Pyrex)
- Centrifuge
- Mixer (Vortex)
- Various pipettes
- HPLC with UV-detector
- Reversed phase C₁₈-column
- Evaporation unit

6. PREPARATION AND STORAGE OF REAGENTS

Reconstitute the **calibrator** in **0.5 ml** reconstitution solution. Take aliquots and store them at -20°C for at least 14 days. Avoid repeated freeze-thaw cycles. The concentration of coenzyme Q10 might have minor changes from lot to lot.

All other test reagents are stable at 2-8 °C, up to the date of expiry stated on the label.

7. PRECAUTIONS

- For in vitro diagnostic use only.
- This product contains human source material which was tested and found to be non-reactive to HBsAg, anti-HIV-1/2, and anti-HCV. Since no method can offer complete assurance that hepatitis B virus, HIV-1/2, HVC or other infectious agents are absent, these reagents should be handled as if potentially infectious.
- Reagents should not be used beyond the expiration date shown on kit label.

8. SPECIMEN COLLECTION AND PREPARATION

EDTA-whole-blood, serum or EDTA-plasma could be used in this assay.

The samples are stable in the dark at 2-8°C for 1 day. For longer storage samples should be frozen at -20°C.

9. ASSAY PROCEDURE

Procedural notes

- The quality control guidelines should be observed.
- Incubation time, incubation temperature and pipetting volumes of the different components are defined by the producer. Any variations of the test procedure, that are not coordinated with the producer, may influence the test results. Immundiagnostik can therefore not be held reliable for any damage resulting from this.
- The assay should always be performed due to the manual which is given in the kit.

Sample and standard preparation

Pipett into 10 ml screw capped glas vials

200 µl sample, standard or control

+

800 µl dilution solution

Mix well (10 s vortex mixer)

Add **1 ml** internal standard

mix well (10 s vortex mixer)

Add **2 ml** extraction solution

Mix for **2 min** on a vortex mixer and centrifuge for **10 min** at 3.000 g. The upper layer is removed and evaporated.

Resuspend the residue in **150 µl** ethanol p.A.

Inject **100 µl** of the resuspension into the HPLC system.

Chromatographic conditions

Column material:	Bischof Prontosil AQ 5 µm
Column dimension:	125 mm x 4 mm
Flow rate:	0.8 - 1.2 ml/min
Detection :	275 nm
Injection volume:	100 µl
Running time:	15 min

Immundiagnostik recommends to use a guard-column to enlarge lifetime of the column.

10. TREATMENT OF THE COLUMN

After analysis the column could be left in the mobile phase. Before use, the system should be equilibrated with ca. 50 ml eluent.


11. RESULTS

Calculation

$$\frac{\text{Peak height patient} \cdot \text{concentration of the calibrator}}{\text{Peak height internal standard patient}} \cdot F = \text{Concentration patient sample}$$

$$F = \frac{\text{Peak height internal standard of the calibrator}}{\text{Peak height of coenzyme Q10 of the calibrator}}$$

Typical chromatogram

12. LIMITATIONS

13. QUALITY CONTROL

Expected values

EDTA-blood: 0.67 – 0.99 $\mu\text{g/ml}$ (mean \pm 2 SD)

EDTA-plasma: 0.83 – 1.43 $\mu\text{g/ml}$ (mean \pm 2 SD)

We recommend, that each laboratory should develop their own normal range. The values mentioned above are only for orientation and can deviate from other published data.

Controls

Control samples or blood pools should be analyzed with each run of calibrators and patient samples. Results generated from the analysis of the control samples should be evaluated for acceptability using appropriate statistical methods. In assays in which one or more of the quality control sample values lie outside the acceptable limits, the results for the patient sample may not be valid.

14. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

Intra-Assay CV:	4.4 % (0.66 µg/ml)	[n = 12]
Inter-Assay VK:	6,6 % (0,31 µg/ml)	[n = 10]
	4,5 % (0,89 µg/ml)	[n = 10]

Linearity

up to 10 µg/ml

Detection limit

0.02 µg/ml

15. DISPOSAL

The mobile phase, extraction solution, internal standard, and ethanol p.A. must be disposed as non-halogenated solvent.

Please refer to the appropriate national guidelines.

16. TROUBLESHOOTING

Problem	Possible reason	Solution
No signal	No or defect connection to evaluation system	Check signalcord and connection
	Detectorlamp is altered	Change lamp
No peaks	Injector is congested	Check Injector
Doublepeaks	Dead volume in fittings and / or column	Renew fittings and / or column
Contaminating peaks	Injector dirty	Clean injector
	Contamination at the head of the column	Change direction of the column and rinse for 30 min at low flow rate (0.2 ml/min) with mobile phase
	Air in the system	Degas pump
	Autosampler vials contaminated	Use new vials or clean them with methanol
Broad peaks, tailing	Precolumn / column exhausted	Use new precolumn / column
Variable retention times	Drift in temperature	Use a column oven
	Pump delivers imprecise	Check pump, degas the system
	System is not in steady state yet	Rinse system mobile phase for 15 min
Baseline is drifting	Detector lamp did not reach working temperature yet	Wait
	Detector lamp is too old	Renew lamp
	System is not in steady state yet	Rinse system mobile phase for 15 min
	Pump delivers imprecise	Check pump, degas the system
Baseline is not smooth	Pump delivers imprecise	Check pump, degas the system
	Detector flowcell is dirty	Clean flow cell

17. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- The test components contain organic solvents. Contact with skin or mucous membranes has to be avoided.
- All reagents in the test package are to be used for research only.
- The reagents should not be used after the date of expiry stated on the label.
- Single components with different lot numbers should not be mixed or exchanged.
- The guidelines for medical laboratories should be observed.
- Incubation time, incubation temperature and pipetting volumes of the different components are defined by the producer. Any variations of the test procedure, that are not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held reliable for any damage resulting from this.