Malondialdehyd

HPLC

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For the determination of Malondialdehyde in plasma and serum

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1. INTENDED USE	15
2. SUMMARY AND EXPLANATION OF THE TEST	15
3. PRINCIPLE OF THE TEST	16
4. MATERIAL SUPPLIED	16
5. MATERIAL REQUIRED BUT NOT SUPPLIED	17
6. PREPARATION AND STORAGE OF REAGENTS	17
7. PRECAUTIONS	17
8. SPECIMEN COLLECTION AND PREPARATION	18
9. ASSAY PROCEDURE	18
PROCEDURAL NOTES	18
SAMPLE AND STANDARD PREPARATION	19
CHROMATOGRAPHIC CONDITIONS	19
10. TREATMENT OF THE COLUMN	20
11. RESULTS	20
CALCULATION	20
TYPICAL CHROMATOGRAM	20
12. LIMITATIONS	21
13. QUALITY CONTROL	21
EXPECTED VALUES	21
CONTROLS	21
14. PERFORMANCE CHARACTERISTICS	22
PRECISION AND REPRODUCIBILITY	22
LINEARITY	22
DETECTION LIMIT	22
15. DISPOSAL	22
16. TROUBLESHOOTING	23
17. REFERENCES	24
18. GENERAL NOTES ON THE TEST AND TEST PROCEDURE	25

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Malondialdehyde

HPLC – Application

*For the determination of Malondialdehyde
in plasma and serum*

Catalogue No.: KC 1900
Package size: 100 determinations
Storage: 2-8 °C standard -20 °C

CE

1. INTENDED USE

The *Immundiagnostik* Assay is intended for the quantitative determination of malondialdehyde in plasma and serum. This Assay is designed for in vitro diagnostic use only.

2. SUMMARY AND EXPLANATION OF THE TEST

In the last years the damage effects of lipid peroxidation products were studied intensively. These are formed when free radicals will overcome the radical-trapping-mechanisms of the body and reacting with unsaturated fatty acids.

The reaction of polyunsaturated fatty acids (PUFA's) with activated oxygen species results in lipidhydroperoxides (primary lipid peroxidation products) which were degraded to secondary lipid peroxidation products like alkanes (e.g. ethane and pentane), aliphatic aldehydes (e.g. malondialdehyde [MDA] and 4-hydroxynonenal [4-HNE]).

Primary and secondary lipid peroxidation products have an influence on a lot of molecules responsible for correct cell function.

So lipidhydroperoxides easily pass the nuclear membrane and can react with nucleic acids. Also proteins can be attacked on their thiol groups changing their stereochemistry and function.

Moreover lipidhydroperoxides interfere with chemical and physical properties of the cell membrane. The fluidity decreases and rigidity increases. The so influenced cell membrane can't maintain their barrier function and intracellular potassium ions leak out as well as intracellular enzymes are lost. If erythrocytes are afflicted haemolysis takes place. In this case haemoglobin can initiate or propagate the lipid peroxidation.

Secondary lipid peroxidation products like MDA or 4-HNE can react with DNA as well, in particular with the bases guanine and adenine. These DNA aberrations lead to erroneous transcriptions and therewith to altered gene products. Peptide bonds are broken through the impact of MDA in proteins. Aldehydes react with amino groups of proteins building Schiff-bases, elementary disturbing correct function of proteins.

All these toxic features of oxidized fatty acids are discussed in the pathogenesis of many diseases and dysfunction of organs. In particular arteriosclerosis, tumor diseases, rheumatic diseases and reperfusion damage of organs after ischemic processes.

3. PRINCIPLE OF THE TEST

The first step in determining malondialdehyde is a sample preparation with derivatisation. The derivatisation reagent transforms malondialdehyde into a fluorescent product. Afterwards the pH is optimized through to addition of a reaction solution. 20µl of the supernatant are injected into the HPLC system.

The separation via HPLC follows an isocratic method at 30°C using a reversed phase column. One run lasts 4 minutes. The chromatograms are recorded by a fluorescence detector. The quantification is performed with the delivered plasma calibrator; the concentration is calculated via integration of the peak heights by the external standard method.

Summary:

Besides many other parameters the advantage of HPLC methode lies in the simultaneous handling of many analytes in a single test. The HPLC system enables even laboratories without experience in high performance liquid chromatography to use this technique for clinical routine determination in a quick and precise manner. Unlike immuno assays with up to six calibrators per test, a one-point calibration is mostly sufficient to calibrate the test system. It is possible to automate the sample application and calculation of the results so that even higher numbers of samples can be handled nearly without control.

4. MATERIAL SUPPLIED

Catalogue No	Kit Components	Quantity
KC 1900lm	Mobile phase	1 x 1000 ml
KC 1900ka	Calibrator, lyophilized	5 vials
KC 1900dl	Derivatisation solution	1 x 100 ml
KC 1900rl	Reaction solution	1 x 50 ml

5. MATERIAL REQUIRED BUT NOT SUPPLIED

- 1.5 ml reaction tubes (Eppendorf)
- Centrifuge
- Various pipettes
- HPLC with Fluorescence-detector
- Reversed phase C₁₈-column - Bischoff ProntoSil Eurobond C18 5 µm, 125 x 4 mm
- Waterbath or heating block for heating at 95 °C

6. PREPARATION AND STORAGE OF REAGENTS

- Reconstitute the **calibrator** in **250 µl** aqua bidest. Take aliquots and store them at -20°C (the reconstituted calibrator is stable for at least 2 weeks at -20°C). Avoid repeated freeze-thaw circles. The concentration of malondialdehyde might have minor changes from lot to lot.
- All other test reagents are stable at 2-8 °C, up to the date of expiry stated on the label.

7. PRECAUTIONS

- For in vitro diagnostic use only.
- This product contain human source material which was tested and found to be non-reactive to HBsAg, anti-HIV-1/2, and anti-HCV. Since no method can offer complete assurance that hepatitis B virus, HIV-1/2, HVC or other infectious agents are absent, these reagents should be handled as if potentially infectious.
- The derivatisation solution contains acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped out immediately with copious quantities of water. Do not breath vapor and avoid inhalation.
- Reagents should not be used beyond the expiration date shown on kit label.

8. SPECIMEN COLLECTION AND PREPARATION

Plasma and serum is suited for this test system. After venipuncture the sample should be stored at 2-8 °C immediately. Samples are stable for at least 24 h at 2-8 °C and 2 weeks at –20 °C.

9. ASSAY PROCEDURE

Procedural notes

- The quality control guidelines should be observed.
- Incubation time, incubation temperature and pipetting volumes of the different components are defined by the producer. Any variations of the test procedure, that are not coordinated with the producer, may influence the test results. Immundiagnostik can therefore not be held reliable for any damage resulting from this.
- The assay should always be performed due to the manual which is given in the kit.

10. TREATMENT OF THE COLUMN

After the analysis the column should be flushed with 15 ml aqua bidest (1 ml/min) and stored in 50% methanol in aqua bidest (ca. 10 ml, flow 0.7 ml/min). Before use, the system should be equilibrated with ca. 15 ml eluent.

Important: The mobile phase must not be recirculated.


11. RESULTS

Calculation

Take into account that the mentioned heights of calibrator and patient are the heights from which the blank value has been subtracted.

$$\frac{\text{peak height patient} \cdot \text{concentration of the calibrator}}{\text{peak height calibrator}} = \text{concentration patient sample}$$

Typical chromatogram

12. LIMITATIONS

Don't use whole blood.

Strong haemolytic and lipaemic samples often show pathological concentrations. We recommend not to measure these samples.

13. QUALITY CONTROL

Expected values

Serum, heparin plasma: $1.97 \pm 0.41 \mu\text{mol/l}$

EDTA plasma: $0.61 \pm 0.24 \mu\text{mol/l}$

We recommend, that each laboratory should develop their own normal range. The values mentioned above are only for orientation and can deviate from other published data.

Controls

Control samples or serum pools should be analyzed with each run of calibrators and patient samples. Results generated from the analysis of the control samples should be evaluated for acceptability using appropriate statistical methods. In assays in which one or more of the quality control sample values lie outside the acceptable limits, the results for the patient sample may not be valid.

14. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

Intra-Assay CV:	6.1 % (0.48 µmol/l)	[n = 12]
	4.8 % (2.06 µmol/l)	[n = 12]
Inter-Assay CV:	6.9 % (0.46 µmol/l)	[n = 12]
	5.7 % (2.13 µmol/l)	[n = 12]

Linearity

up to 100 µmol/l

Detection limit

0.15 µmol/l

15. DISPOSAL

The mobile phase, derivatisation solution and reaction solution must be disposed as non-halogenated solvent.

Please refer to the appropriate national guidelines.

16. TROUBLESHOOTING

Problem	Possible reason	Solution
No signal	No or defect connection to evaluation system	Check signalcord and connection
	Detectorlamp is altered	Change lamp
No peaks	Injector is congested	Check Injector
Doublepeaks	Dead volume in fittings and / or column	Renew fittings and / or column
Contaminating peaks	Injector dirty	Clean injector
	Contamination at the head of the column	Change direction of the column and rinse for 30 min at low flow rate (0.2 ml/min) with mobile phase
	Air in the system	Degas pump
	Autosampler vials contaminated	Use new vials or clean them with methanol
Broad peaks, tailing	Precolumn / column exhausted	Use new precolumn / column
Variable retention times	Drift in temperature	Use a column oven
	Pump delivers imprecise	Check pump, degas the system
	System is not in steady state yet	Rinse system mobile phase for 15 min
Baseline is drifting	Detector lamp did not reach working temperature yet	Wait
	Detector lamp is too old	Renew lamp
	System is not in steady state yet	Rinse system mobile phase for 15 min
	Pump delivers imprecise	Check pump, degas the system
Baseline is not smooth	Pump delivers imprecise	Check pump, degas the system
	Detector flowcell is dirty	Clean flow cell

17. REFERENCES

Draper H.H., Hadley M. (1990). A review of recent studies on the metabolism of exogenous and endogenous malondialdehyde. *Xenobiotika* 20; 9; 901-907.

Griesmacher A. et al. (1995). Enhanced serum levels of thiobarbitric-acid-reactive substances in diabetes mellitus. *Am J Med* 98; 469-475.

Valenzuela A (1990). The biological significance of malondialdehyde determination in the assessment of tissue oxidative stress. *Life sciences* 48; 301-309

18. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- The test components contain organic solvents. Contact with skin or mucous membranes has to be avoided.
- All reagents in the test package are to be used for research only.
- The reagents should not be used after the date of expiry stated on the label.
- Single components with different lot numbers should not be mixed or exchanged.
- The guidelines for medical laboratories should be observed.
- Incubation time, incubation temperature and pipetting volumes of the different components are defined by the producer. Any variations of the test procedure, that are not coordinated with the producer, may influence the results of the test. Immundiagnostik AG can therefore not be held reliable for any damage resulting from this.