

BÜHLMANN Allergen List



INSECT VENOMS	
BAG2-I1	HONEY BEE VENOM
BAG2-I1CHK	HONEY BEE VENOM CHECK <i>new</i>
BAG2-I3	YELLOW JACKET VENOM
BAG2-I3CHK	YELLOW JACKET VENOM CHECK <i>new</i>
BAG2-I4	PAPER WASP VENOM
BAG2-I75	EUROPEAN HORNET VENOM
DRUGS	
Antibiotics	
BAG2-C1	PENICILLIN G
BAG2-C11	BENZYLPENICILLOYL-POLYLYSIN, PPL
BAG2-C12	BENZYLPENICILLIN +-ACID, MDM
BAG2-C2	PENICILLIN V
BAG2-C3	CEPHALOSPORIN C
BAG2-C31	CEFAMANDOLE
BAG2-C32	CEFAZOLIN
BAG2-C33	CEFUROXIME
BAG2-C34	CEFACTOR <i>new</i>
BAG2-C35	CETRIAXONE <i>new</i>
BAG2-C36	LEVOFLOXACIN <i>new</i>
BAG2-C61	SULFAMETHOXAZOLE
BAG2-C62	TRIMETHOPRIM
BAG2-C75	TETRACYCLINE
BAG2-C81	CIPROFLOXACIN
BAG2-C203	AMPICILLIN (EK-CAST)
BAG2-C203F	AMPICILLIN (FK-BAT)
BAG2-C204	AMOXICILLIN
BAG2-CCLA	CLARYTHROMYCIN
BAG2-CRIF	RIFAMPICIN
Analgesics	
BAG2-C51	LYS-ASPIRIN
BAG2-C52	DICLOFENAC
BAG2-C53	IBUPROFEN
BAG2-C54	INDOMETHACIN
BAG2-C55	ACETAMINOPHEN
BAG2-C56	MEFENAMIC ACID
BAG2-C57	PHENYL BUTAZONE
BAG2-C58	PROPYPHENAZONE
BAG2-C59	DIPYRONE / METAMIZOLE
BAG2-CNAP	NAPROXEN
Myorelaxants	
BAG2-CATR	ATRACURIUM
BAG2-CCAT	CISATRACURIUM <i>new</i>
BAG2-CMIV	MIVACURIUM
BAG2-CPAN	PANCURONIUM
BAG2-CPRO	PROPOFOL
BAG2-CROC	ROCURONIUM
BAG2-CSUX	SUXAMETHONIUM
BAG2-CVEC	VECURONIUM
Local Anesthetics	
BAG2-CART	ARTICAINE <i>new</i>
BAG2-CBUP	BUPIVACAINE <i>new</i>
BAG2-CMEP	MEPIVACAINE <i>new</i>
BAG2-CLID	LIDOCAINE
Contrast Media	
BAG2-CIOBI	IOBITRIDOL <i>new</i>
BAG2-CIODI	IODIXANOL <i>new</i>
BAG2-CIOHE	IOHEXOL <i>new</i>
BAG2-CIOME	IOMEPROL <i>new</i>
BAG2-CIOPA	IOPAMIDOL <i>new</i>
BAG2-CIOPR	IOPROMIDE <i>new</i>
BAG2-CIOXA	IOXAGLATE <i>new</i>
BAG2-CIOXI	IOXITALAMATE <i>new</i>

Antiseptics	
BAG2-CCHX	CHLORHEXIDINE
ENVIRONMENTAL/OCCUPATIONAL	
BAG-K79	PHTALIC ANHYDRIDE/PHTALIC ACID
BAG-K80	FORMALDEHYDE
BAG-K82	LATEX
BAG-K85	CHLORAMINE-T
BAG-K87	ALPHA-AMYLASE
INHALANTS	
Grasses	
BAG-G2	BERMUDA GRASS
BAG-G3	ORCHARD GRASS
BAG-G5	PERENNIAL RYE GRASS
BAG-G6	TIMOTHY GRASS
BAG-G12	CULTIVATED RYE GRASS
WEEDS	
BAG-W1	COMMON RAGWEED
BAG-W6	MUGWORT
BAG-W9	RIBWORT/PLANTAIN
BAG-W19	PARIETARIA OFFICINALIS
Trees	
BAG-T3	COMMON BIRCH
BAG-T4	HAZEL
BAG-T7	OAK
BAG-T9	OLIVE
Molds	
BAG-M1	PENICILLIUM NOTATUM
BAG-M2	CLADOSPORIUM HERBARUM
BAG-M3	ASPERGILLUS FUMIGATUS
BAG-M5	CANDIDA ALBICANS
BAG-M6	ALTERNARIA TENUIS
Mites/Animals/Insects	
BAG-D1	D. PTERONYSSINUS
BAG-D2	D. FARINAE
BAG-D70	ACARUS SIRO
BAG-E1	CAT EPITHELIUM
BAG-E2	DOG EPITHELIUM
BAG2-CUL	CULICOIDES ANTIGEN (RUO*, HORSE APPLICATION)
FOOD	
Egg	
BAG-F1	EGG WHITE
BAG-F75	EGG YOLK
Milk/Milk Products	
BAG-F2	MILK (PS75, DEFATTED)
BAG-F76	ALPHA-LACTALBUMIN
BAG-F77	BETA-LACTOGLOBULIN
BAG-F78	CASEIN
Fish & Meat	
BAG-F3	CODFISH
BAG-F23	CRAB
BAG-F24	SHRIMP
BAG-F258	SQUID
BAG-F290	OYSTER
BAG-P4	ANISAKIS
BAG-F26	PORK
BAG-F27	BEEF
Fruits/Vegetables	
BAG-F25	TOMATO
BAG-F31	CARROT
BAG-F33	ORANGE
BAG-F95	PEACH
BAG-F85	CELERY

Allergens/Allergen Mixes



Seeds/Beans/Nuts

BAG-F10	SESAME
BAG-F13	PEANUT
BAG-F14	SOYBEAN
BAG-F17	HAZELNUT
BAG-F20	ALMOND
BAG-F202	CASHEW NUT

Cereales

BAG-F4	WHEAT
BAG-F5	RYE FLOUR
BAG-F6	BARLEY FLOUR
BAG-F7	OAT FLOUR
BAG-F8	MAIZE
BAG-F81	MAIZE STARCH <i>new</i>
BAG-F9	RICE
BAG-F45	BAKER'S YEAST

FOOD ADDITIVES

BAG2-C103	TARTRAZINE
BAG2-C105	METHYLENE BLUE <i>new</i>
BAG2-C111	SODIUM-BENZOATE
BAG2-C112	SODIUM-NITRITE
BAG2-C113	POTASSIUM-METABISULFITE
BAG2-C114	SODIUM-SALICYLATE
BAG2-CE104	QUINOLINE YELLOW
BAG2-CE110	SUNSET YELLOW FCF
BAG2-CE122	CHROMOTROPE B
BAG2-CE123	AMARANTH
BAG2-CE124	NEW COCCINE
BAG2-CE127	ERYTHROSINE
BAG2-CE131	PATENT BLUE V
BAG2-CE132	INDIGO CARMINE
BAG2-CE151	BRIILLIANT BLACK BN
BAG2-CE466	CARBOXYMETHYLCELLULOSE
BAG2-CE621	GLUTAMATE

RECOMBINANT ALLERGENS*

BAG2-ZEAM12	ZEAM 12 (MAIZE PROFILIN)
BAG2-TRIA12	TRIA 12 (WHEAT PROFILIN)
BAG2-FELD1	FEL D 1 (CAT MAJOR ALLERGEN)
BAG2-MALD1	MAL D 1 (APPLE MAJOR ALLERGEN)

*(RUO, Research use only)

ALLERGEN MIXES

BAG-DX1	STORAGE MITE
ACARUS SIRO (D70), LEPIDOGLYPHUS DESTRUCTOR, TYROPHAGUS PUTRESCENTIA, GLYCYPHAGUS DOMESTICUS	

BAG-GX1	6-GRASS
ORCHARD GRASS (G3), MEADOW FESCUE GRASS, PERENNIAL RYE GRASS (G5), TIMOTHY GRASS (G6), MEADOW GRASS, VELVET GRASS	

BAG-WX1	RAGWEED
COMMON RAGWEED (W1), GIANT RAGWEED	

BAG2-C101	FOOD COLORANTS I
QUINOLINE YELLOW (CE104), SUNSET YELLOW FCF (CE110), CHROMOTROPE B (CE122), AMARANTH (CE123), NEW COCCINE (CE124)	

BAG2-C102	FOOD COLORANTS II
ERYTHROSINE (CE127), PATENT BLUE V (CE131), INDIGO CARMINE (CE132), BRILLIANT BLACK BN (CE151)	

BAG-FX1	FOOD MIX
EGG WHITE (F1), EGG YOLK (F75), MILK (F2), CODFISH (F3), WHEAT (F4), PEANUT (F13), SOYBEAN (F14), HAZELNUT (F17), SHRIMP (F24)	

BAG-IX1	INHALANT MIX
6-GRASS MIX (GX1), CULTIVATED RYE GRASS (G12), COMMON BIRCH (T3), HAZEL (T4), MUGWORT (W6), RIBWORT/PLANTAIN (W9), ALTERNARIA TENUIIS (M6), D. PTERONISSINUS (D1), D. FARINAE (D2), CAT EPITHELIUM (E1), DOG EPITHELIUM (E2)	

BAG2-CNBX	NEUROMUSCULAR BLOCKERS MIX <i>new</i>
ATRACURIUM (CATR), CISATRACURIUM (CCAT), MIVACURIUM (CMIV), PANCURONIUM (CPAN), PROPOFOL (CPRO), ROCURONIUM (CROC), SUXAMETHONIUM (CSUX), VECURONIUM (CVEC)	

Resuspend allergens with 250 µl Stimulation Buffer (STB)

Up to 4 stimulations per allergen vial

Store unopened allergens at 2-8°C until exp. date

Use a fresh allergen vial each time a new stimulation is performed; except Bee & Wasp allergens:

Reconstituted Bee & Wasp allergens (BAG2-I1 & BAG2-I3) are stable for 1 month stored at -20°C

For further information refer to the Allergen Data Sheets available on the BÜHLMANN webpage www.buhmannlabs.ch. Search term: Allergens



BÜHLMANN Laboratories AG
Baselstrasse 55
CH-4124 Schönenbuch/Basel
Switzerland

Phone +41 61 487 12 12
Fax orders +41 61 487 12 99
info@buhmannlabs.ch
www.buhmannlabs.ch



Ordering code:	
EK-CAST	192 wells
EK-CAST5	480 wells
FK-CCR	100 tests
FK-BAT	100 tests
FK-BAT2	200 tests

CAST® is a registered trademark of BÜHLMANN Laboratories AG