

Now Available!



Vitrification

freeze your future™

Irvine Scientific's Vitrification System is a state-of-the-ART cryopreservation method that can be used for oocytes and embryos.

Vitrification offers many benefits compared to conventional freeze/thaw techniques:

- Simple, straight forward procedure
- Our unique CryoTip™ is the only closed system available
- The entire procedure takes just a few minutes, not hours
- No special equipment is required
- Greater than 90% recovery post-thaw, yielding greater than 50% pregnancy rates*

CryoTips™ are specially designed, drawn plastic straws with ultra-fine tips and protective metal cover sleeves. CryoTips™ provide a closed system device for Vitrification. This ensures that the specimens do not have direct contact with liquid nitrogen or other potential sources of contamination.

Vit Kit™ - Freeze

- 2 vials of Equilibration Solution
- 4 vials of Vitrification Solution
- 10 CryoTips™
- 2 Connectors
- Demonstration video

Vit Kit™ - Thaw

- 2 vials of Thawing Solution
- 2 vials of Dilution Solution
- 2 vials of Washing Solution
- 2 Connectors
- Demonstration video

*Masaashige Kuwayama et al. *RBM Online*, Volume 11 - Accepted articles unassigned to issue.
Abstract Accepted April 2005.

Irvine Scientific

2511 Daimler Street, Santa Ana, CA 92705

Customer Service: 1 800 577 6097 • Fax: 1 949 261 6522

Telephone: 1 949 261 7800 • Technical Service: 1 800 437 5706

www.irvinesci.com • vitkit_info@irvinesci.com

European Representative
MPIE
Schuurweg 13a
5145 NP Waalwijk, The Netherlands
CE 0050 See instructions
for use.



IrvineScientific®

Grow With Us